

## FIBROBLAST MEDIUM-PLUS

---

<b>Product Type:</b>	Fibroblast Medium-Plus
<b>Catalog Number:</b>	P60166

---

### **Product Description**

Fibroblast Medium-Plus (FM-Plus) is a complete medium designed for optimal growth of normal human cardiac and kidney fibroblasts in vitro. It is a sterile, liquid medium which contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals and a low concentration of fetal bovine serum (5%). The medium is HEPES and bicarbonate buffered and has a pH of 7.4 when equilibrated in an incubator with an atmosphere of 5% CO<sub>2</sub>/95% air. The medium is formulated (quantitatively and qualitatively) to provide a defined and optimally balanced nutritional environment that selectively promotes proliferation and growth of normal human cardiac fibroblasts in vitro.

### **Components**

- 500 ml of Basal Medium
- 25 ml of Fetal Bovine Serum (FBS)
- 5 ml of Fibroblast Growth Supplement 2 (FGS-2)
- 5 ml of penicillin/streptomycin solution (P/S solution)

### **Prepare for use**

Thaw FGS-2, FBS and P/S solution at 37°C. Gently tilt the FGS tube several times during thawing to help the contents dissolve. **Make sure the contents of the supplement are completely dissolved into solution before adding to the medium.** Rinse the bottle and tubes with 70% ethanol, and then wipe to remove excess. Remove the cap, being careful not to touch the interior threads with fingers. Add FGS, FBS and P/S solution into basal medium in a sterile field, mix well and then the reconstituted medium is ready for use. Since several components of this medium are light-labile, it is recommended that the medium not be exposed to light for lengthy periods of time. If the medium is warmed prior to use, do not exceed 37°C. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

### **Caution**

*If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.*