

# **Experimental Protocol**

# FluoHitSeeker FSHR

Catalog #: P30117-F

#### 1. Introduction

FluoHiTSeeker reporter assays are used widely to investigate cellular signaling pathways and as high-throughput screening tools for drug discovery. The FluoHiTSeeker FSHR Cell Line is a clonal derivative of Human Embryonic Kidney 293 (HEK 293) cells. These cells contain a turboGFP gene (tGFP) under the control of a minimal promoter with cAMP Response Element (CRE). CRE is the DNA-binding sequence for the transcription factor CRE binding protein (CREB), which is responsible for the regulation of a variety of biological functions including cell proliferation, circadian rhythms and memory. FluoHiTSeeker FSHR Cell Line is designed for High throughput screening (HTS) analysis of receptor response that results in modulation of CREB activities. Elevation of the intracellular cAMP level activates cAMP response element binding protein (CREB) to bind CRE and induces the expression of the fluorescent protein turboGFP.

#### 2. Product Components and Storage Conditions

Product: FSHR FluoHiTSeeker

cat.nº. P30117-F

Size: 2 vial 3x106cells in Freezing Media

## 3. Biological Activity:

This cell line is validated for cellular response to stimulation by Follicle-stimulating hormone.

#### Warranty

## Cell Line Stability

Cells may undergo genotypic changes resulting in reduced responsiveness over time in normal cell culture conditions. Genetic instability is a biological phenomenon that occurs in all stably transfected cells. Therefore, it is critical to prepare an adequate number of frozen stocks at early passages.

#### Mycoplasma testing

The cell line has been screened using the PCR-based Venor™GeM Mycoplasma Detection kit (Sigma-Aldrich) to confirm the absence of Mycoplasma species.

#### Storage

Immediately upon receipt, store in liquid nitrogen.



## 4. Materials to Be Supplied by the User

Recommended Reagents
DMEM, high glucose (Sigma-Aldrich D6429)
Fetal bovine serum (FBS)
DPBS with calcium and magnesium (Sigma Aldrich D8662)
Opti-MEM (Life technologies 31985-070)

## 5. Supplies and Equipment

96-well assay plate
Tissue culture flasks
Class II biological safety cabinet
Hemacytometer
Incubator humidified 37°C, 5% CO2
Inverted microscope
Fluorimetre

## 6. Experimental protocol

- 1. Harvest FSHR FluoHiTSeeker-HEK293 cells from culture in growth medium and seed cells at a density of ~50,000 cells per well in 200 μl of culture medium into a white clear-bottom 96-well microplate.
- 2. Incubate cells at 37°C in a CO<sub>2</sub> incubator overnight (~ 16-18 hours).
- 3. Remove the media and add threefold serial dilution of Follicle-Stimulating Hormone in Opti-MEM to stimulated wells. Add Opti-MEM including the vehicle of the compounds to the unstimulated control wells. Set up each treatment in at least triplicate.
- 4. Incubate the plate at 37°C in a CO2 incubator for 16-24 hours.
- 5. Remove the Assay medium and replace by 100 µl of DPBS with calcium and magnesium. Read the plate using the appropriate filter for the tGFP protein fluorescent signal (excitation/ emission max = 482/502 nm).
- Data Analysis: Subtract average background of fluorescence (FSH-free control wells) from fluorescence reading of all wells.